

Observing Life As It Happens



**QUICKSTART GUIDE:**  
**WIDEFIELD WF3**  
**Zeiss Cell Observer**  
**Live Cell Imaging System**  
(SAF, ROOM 409)



|   |   |
|---|---|
| Startup procedure .....                         | 2 |
| Viewing your sample .....                       | 2 |
| Adjusting brightfield Köhler illumination ..... | 3 |
| Setting user-defined pre-sets .....             | 3 |
| Image Acquisition .....                         | 3 |
| Saving .....                                    | 4 |
| Additional Methods .....                        | 4 |
| Software Autofocus .....                        | 4 |
| Z stack .....                                   | 5 |
| Tile and Positions .....                        | 5 |
| Time series .....                               | 5 |
| Heating & CO <sub>2</sub> units .....           | 6 |
| Shutdown procedure .....                        | 6 |

# Startup procedure

- All switches on wall to on
- Check camera is on (left hand side)
- Start PC and login
- Run Zen blue and click on “Zen Pro”
- If you are the first user after switching on – the stage will want to calibrate. Select the 5x objective on the TFT controller and calibrate



## TFT touchscreen control


Can be used to:

- Select objective
- Switch on/off Transmitted light
- Switch on/off Reflected light
- Select reflector cube
- Control heating and CO2
- Additional focussing control
- Other microscope controls




## Viewing your sample

- In the software select **locate** tab
- Select **BF** (brightfield)
- Turn on Transmitted light either in software or using the TFT touchscreen



## Adjusting brightfield Köhler illumination


- Select 10x objective either with the TFT touch screen or software
- you can move your specimen in xy-direction with the joystick, holding the button on top of the joystick makes the stage moving faster
- focus on your sample
- adjust brightness with large black wheel on front / bottom of microscope)
- make sure condenser iris is not completely closed (two buttons on left top of condenser turret)
- completely close field iris (top black wheel) (if the image turns completely black, reopen until you see some light, proceed to next step and close again after you have focussed the iris)
- focus iris (black condenser focus wheel)
- centre iris (2 silver Köhler screws)
- reopen field iris until the edge is just not visible anymore




## Image Acquisition

- Select the Acquisition tab (make sure “show all tools” checkbox is ticked)
- Either open and image and “Re-Use” or load one of the saved experiment setups (just below the Acquisition tab – Experiment Manager)
  - WF3 Basic 5 channel (Dapi, GFP,TRITC/Texas Red, CY5, Brightfield)
  - WF3 Quad (Fast acquisition - Brightfield, Dapi, GFP,TRITC/Texas Red, CY5)
  - WF3 Colour Camera
  - WF3 Flash BF (Flash camera brightfield only – live cell imaging)
  - WF3 CFP/FRET/YFP (request filter set)
- To adjust Channels use the “Channels” tab
- Click “Live”
- Select each Channel in turn – adjusting LED power (using the Lumencor button) and exposure time (range indicator button is below image window)
- Additional camera setting are in Camera tab
- To capture image click Snap


Select Channel



Change Camera exposure time



- Settings can be saved for future use




Select New

Type in name

Select Save

# Lumencor LED's

Click the button above the image window and this will open the LED intensity controller




These settings are not stored with the image so need to be remembered (use the windows snipping tool to capture and save the settings)

## Saving

- Images are stored in your temporary folder until you “save” them
- To save the image – select it and press the save icon.
- Select location and name

## Additional Methods

Select check box for Z stack, Tile and Time Series. This will activate the respective tabs  
To capture images using these methods use the Start Experiment button



## Software Autofocus


Select Focus Strategy tab  
Options will depend on what other functions are selected eg Z stack and Tile

- Absolute fixed Z position
- Software autofocus
- Definite focus

- Software autofocus as a reference for Definite focus
- Definite focus as a reference for Software autofocus
- Local focus
- Local focus for software autofocus


## Z stack

- Set First and Last positions
- Set interval
- Start experiment



## Tile and Positions

- Setup tile or positions
- Start experiment
- (ask for more instruction)



## Time series

- Setup duration and cycles
- Start experiment




## For very fast acquisition single channel

- Load WF3 Flash BF protocol
  - Select fluorochrome
  - Select "Speed"
- Set the camera settings as required
- Go "Live" and set exposure and focus
- Setup duration under "Time Series" and tick "use Camera streaming if possible"
- Run "Start experiment"

## Heating & CO<sub>2</sub> units

- Found to the right of the microscope
- Switch on heating at least 1 hour before required
- For CO<sub>2</sub> open cylinder in room 408
- Unclip hoseclips and adjust flow



## Shutdown procedure

- check if anyone is booked after you within 2 hours

### If nobody is booked within two hours:

- remove specimen
- clean oil objective lenses with lens tissue
- Important: close incubation chamber completely (→ protects from dust)
- update time in Sharepoint
- secure your data (e.g. copy them to the server)
- shut down hardware in reverse order from startup

### If someone's booked within two hours:

- update usage in Sharepoint
- remove your samples
- clean objective lenses with fresh lens tissue and close incubation chamber
- clear up the desk
- secure your data (e.g. copy them to the server)
- log off